

Sleep Therapy I Ventilation I Patient Interface

## Sleep Therapy

Sleep therapy solutions for every patient



# Sleep Therapy

Sleep therapy solutions for every patient



## Tailor-made therapy for every patient

Doctors working in many different medical specialties are diagnosing Sleep-Disordered Breathing (SDB) with increasing frequency. The consequences of SDB span the range of daytime sleepiness with decreased productivity and high blood pressure to myocardial infarction and stroke. Patients with heart failure, for example, face a higher risk of suffering from SDB. For decades we have

been developing reliable, specialized technologies for use in sleep medicine. Our focus is on device safety, lasting therapy success and ease of use. The result is products that give the patient effective treatment and, through intuitive operation, product compatibility and easy maintenance, help you to simplify your work processes and allow you to concentrate on the patient.

FOT / FBT To distinguish central from obstructive apnea

SCOPES LAB / ST /

## **Indications**

- Central sleep apnea of Cheyne-Stokes variety associated with heart failure
- and with central, mixed or complex sleep apnea

## Type of therapy

## **AcSV** therapy

Ventilation by means of AcSV combined with automatic pressure adjustment to patient's needs during current breath; automatic or pre-set backup frequency

## Therapy devices

## prismaCR

The sleep medicine therapy solution for complex SDB (Sleep-Disordered Breathing)



- Obstructive mixed or complex sleep apnea and
  - High and/or changing pressure needs
  - Poor CPAP compliance
  - Central apnea
- Nocturnal hypoventilation (as with OHS)
- Respiratory insufficiency, e.g., coprevalent COPD (overlap)

## **BILevel ST therapy**

Ventilation with varied pressure levels during inspiration and during early and late expiration phases; depending on the therapy goal (three SCOPES for simple and fast setting). Backup frequency for additional safety.

## prisma30ST

BILevel ST therapy device with adjustable target volume



### prisma25ST

BILevel ST device with innovative pre-settings concept (SCOPES) and optional automatic therapy monitoring



- Obstructive mixed or complex sleep apnea and
  - High and/or changing pressure needs
  - Poor CPAP compliance

## **BILevel S therapy**

Ventilation with varied pressure levels during inspiration and during early and late expiration phases

## prisma25S

BILevel S device combined with the advantages of APAP



- Obstruction of the upper airways
- Fixed or variable pressure needs depending on position or sleep stage
- Poor CPAP compliance
- OSA with varying pressure needs

## **APAP** therapy

Ventilation with constant pressure that automatically adjusts to the needs of the patient to provide greatest therapy flexibility. As an option, CPAP and APAP can be combined with the intelligent pressure relief softPAP.

## prisma20A

Therapy flexibility for OSA (Obstructive Sleep Apnea)



## SOMNObalance e

For balanced APAP therapy



- Obstruction of the upper airways
- Fixed pressure needs

## **CPAP** therapy

Ventilation with constant pressure

## SOMNOsoft 2e

Comfortable and effective CPAP therapy



**OPP** (Obstructive Pressure Peak) to distinguish central from obstructive apnea



## Overview of sleep therapy devices

## prismaLINE - A platform to trust

## prismaLINE is the intelligent system solution for treatment of respiratory disorders.

The prismaLINE represents a new dimension in operating convenience. The combination of state-of-the-art touch screen and smart Graphical User Interface makes possible unprecedented intuitive and therefore fast usage.

		0 42	0 0 0	06.42	
		prismaLAB	prismaCR	prisma30ST	
Modes	CPAP	<b>V</b>	V	<b>V</b>	
	APAP (autoCPAP)	<b>V</b>	-	<b>V</b>	
	BILevel S, auto S	<b>V</b>	-	<b>V</b>	
	BILevel S/T, auto S/T	<b>V</b>	_	<b>V</b>	
	BILevel T	<b>V</b>	_	<b>V</b>	
	AcSV*	<b>V</b>	V	_	
	aPCV	<b>V</b>	-	<b>V</b>	
Stabilization	Obstruction detection FOT/FBT	<b>V</b>	<b>V</b>	<b>V</b>	
	Adaptive pressure regulation to eliminate obstructions (autoEEPAP, APAP)	<b>V</b>	<b>V</b>	<b>V</b>	
	Pressure support	adaptive, fixed	adaptive	fixed	
	Optional autoF	<b>V</b>	<b>V</b>	<b>V</b>	
Synchron- ization	Adjust target volume	<b>V</b>	_	<b>V</b>	
	Automatic trigger technology	<b>V</b>	<b>V</b>	<b>V</b>	
	Automatic pressure rise	Ramp IN, Ramp EX	auto	Ramp IN, Ramp EX	
Compliance and Comfort	softSTART	<b>V</b>	<b>V</b>	<b>V</b>	
	Mask test	<b>V</b>	<b>V</b>	<b>V</b>	
	Leak display during operation	<b>V</b>	<b>V</b>	<b>V</b>	
	TRILevel option	<b>V</b>	<b>V</b>	<b>V</b>	
	Pressure relief softPAP	in CPAP/APAP modes	in CPAP mode CPAP	in CPAP/APAP modes	
	autoSTART-STOP	<b>V</b>	<b>V</b>	<b>V</b>	
Therapy Info for Experts	Signals for PSG feed	8 signals	8 signals	8 signals	

This exceptional operating concept is found throughout the prismaLINE portfolio. Every prismaLINE device provides the user with ideal support for all kinds of work processes and for connectivity and accessories too. prismaLINE offers the

right device for every type of therapy in Weinmann quality: precise, reliable and unheard-of stillness. prismaLINE: Synonymous with modern sleep therapy Made in Germany.

0 0	0	0 No. 12				
prisma25ST	prisma25S	prisma20A	SOMNO <i>balanc</i> e e	SOMNOsoft 2e		
<b>/</b>	<b>V</b>	<b>V</b>	<b>V</b>	V		
<b>V</b>	V	<b>V</b>	<b>V</b>	_		
<b>V</b>	<b>V</b>	_	_	_		
<b>V</b>	_	-	_	_		
<b>V</b>	_	-	_	_		
-	_	-	_	_		
-	_	-	-	_		
<b>V</b>	<b>V</b>	<b>v</b>	OPP	OPP		
<b>V</b>	<b>V</b>	<b>v</b>	<b>V</b>	_		
fixed	fixed	_	_	_		
<b>V</b>	_	-	-	_		
_	_	_	_	_		
<b>V</b>	<b>V</b>	_	-	_		
<b>✓</b> Ramp IN	<b>V</b>	_	_	_		
<b>V</b>	<b>v</b>	<b>v</b>	<b>V</b>	<b>V</b>		
<b>V</b>	<b>V</b>	<b>v</b>	<b>V</b>	V		
<b>V</b>	<b>V</b>	<b>v</b>	<b>V</b>	<b>V</b>		
<b>V</b>	<b>V</b>	-	_	_		
n CPAP/APAP modes	in CPAP/APAP modes	V	<b>V</b>	V		
<b>V</b>	<b>V</b>	<b>v</b>	<b>V</b>	V		
8 signals	8 signals	8 signals	2 from 5 signals	2 from 5 signals		
* AcSV=Anticyclic Servoventilation						

<sup>\*</sup> AcSV=Anticyclic Servoventilation



## prismaLINE accessories

## prismaLINE meets all needs with numerous add-ons

All prismaLINE models look identical on the outside and offer standardized expansion interfaces. To modify devices for specific uses, the user can connect different modules to any device and combine devices with each other. With so much flexibility, doctors and patients can be sure that the ideal therapy solution will be found, no matter how complex the clinical picture is.



## Compatible and futuristic

Data transmission can be made in any number of ways, all of which are based on proven industry standards that are manufacturer-independent. The therapy devices are equipped with interfaces for SD cards and USB for this purpose. Expansion modules for LAN and PSG are available; devices can be refitted for telemedicine via GSM.\*





prismaAQUA Humidifier WM 29680



prismaCONNECT Connection module WM 29670



prismaPSG PSG module WM 29690

## Micro USB connection cable WM 35160

## Converter 24 V, Lorry, WM 24617 12 V, Auto, WM 24616



**Breathing tube** Ø 15 mm, WM 29988 Ø 22 mm, WM 24445 Ø 22 mm, autoclavable, WM 24667



e.g., JOYCEone Ventilation mask WM 25260 (nasal) WM 25290 (full face)

**Bacteria filter** WM 24476

## prismaHYBERNITE Heated breathing tube Ø 22 mm, WM 29067

## Set 12 Pollen filter WM 29652

Power failure detector WM 29650



## prismaLABset-PSG

With prismaLAB device, prismaAQUA, prismaBAG butler (WM 29978), prismaCONNECT, prismaTS software with USB data cable and prismaPSG

Complete titration set WM 29996



## prismaTS and prismaTSlab

## The accompanying software solutions

From a quick summary to detailed therapy analysis, prismaTS and prismaTSlab provide relevant data for every work process.

Therapy software complete with USB data cable, WM 93335

<sup>\*</sup>Legal regulations limit availability to certain countries.

## SOMNO devices accessories

## Compatible, efficient SOMNO accessories

Accessories include the external leak alarm SOMNOalarm, perfect-fit ventilation masks like JOYCEone and the heated tube system HYBERNITE Superday. With compatible accessories that are tailored exactly to the needs of patients and medical staff, treatment with SOMNOsoft 2e and SOMNObalance e is even more effective and successful.



JIVINO accessories



**SOMNO***aqua* Humidifier WM 24403



SOMNOalarm External leak alarm WM 29750



O<sub>2</sub> safety valve Safe oxygen feed WM 24042



HYBERNITE Superday Heated tube system Ø 22 mm, WM 29045



e.g., JOYCEone Ventilation mask WM 25260 (nasal) WM 25290 (full face)



Breathing tube Ø 22 mm, WM 24445 Ø 22 mm, autoclavable, WM 24667

Bacteria filter WM 24476



WEINMANNsupport PC software for settings and analysis WM 93301

**DC Adapter** 12 V / 24 V WM 24469

**USB Converter cable** WM 93321

 $You'll\ find\ additional\ information\ about\ our\ the rapy\ solutions,\ accessories\ and\ mask\ systems\ at\ \textbf{weinmann-medical.com}$ 

Weinmann Geräte für Medizin GmbH + Co. KG develops, produces and markets diagnostic and therapeutic solutions for sleep medicine and ventilation.

Products from the Hamburg-based, family-owned company are tailored exactly to the needs of patients and experts. Before market launch, every product goes through a certified validation process in research and development, procurement, production and quality management.

Weinmann, a Löwenstein Group company since 2013, is among the Top 3 in many of its demanding international markets, thanks to ideas, experience and innovation. Thanks too to the decades of trust that Weinmann's partners around the world have placed in our proven products Made in Germany.

weinmann-medical.com

### Germany

Weinmann Geräte für Medizin GmbH + Co. KG Postfach 540268 ■ 22502 Hamburg Kronsaalsweg 40 ■ 22525 Hamburg

E: info@weinmann-medical.com weinmann-medical.com

T: +49 (0)40 547 02-0 Reception

F: +49 (0)40 547 02-461 Reception

T: +49 (0)40 547 02-100 Customer Service

Center for Production, Logistics, Service Weinmann Geräte für Medizin GmbH + Co. KG Siebenstücken 14 ■ 24558 Henstedt-Ulzburg

## France

Weinmann S.A.S. – Paris-Igny T: +33 (0)169 355320 weinmann-france.fr



## Russia

Weinmann SPb GmbH – Moscow T: +7 (812)633 3082 weinmann.ru



## China

Weinmann GmbH + Co. KG (Shanghai Rep. Office) T: +86 21 32 26 26 98 weinmann.cn



